

Product Datasheet

**Recombinant SARS-CoV-2 (COVID-19)
Nucleocapsid ScFv-Fc Antibody [4A6]
(orb1274613)**

Description

Recombinant SARS-CoV-2 (COVID-19)


Species/Host

Human

Reactivity

Virus

Conjugation

Unconjugated

Tested Applications

ELISA

Immunogen

full length recombinant Nucleocapsid Protein

Target

N

Preservatives

PBS, pH 7.4

Form/Appearance

Liquid

Concentration

5.0 mg/ml (+/- 10 %)

Storage

 24 months from manufacturing at -20 °C;
1 month, 2 to 8 °C under sterile conditions.

Note

For research use only

Application notes

Immunoreactivity: 80-120 % compared to the reference sample in an ELISA test.

Isotype

Human IgG1

Clonality


Recombinant

Uniprot ID
[P0DTC9](#)
NCBI
[QHD43423](#)
Dilution Range


Immunoreactivity: 80-120 % compared to the reference sample in an ELISA test.

Expiration Date

12 months from date of receipt.

 Immobilized
2019-nCoV
Nucleocapsid
Prote...


Stability is confirmed by binding abilit...



Please note that pair recommendations ar...