


Product Datasheet

SARS-CoV-2 Spike NTD Antibody [2146] (orb1274626)

Description	SARS-CoV-2 Spike NTD Antibody [2146]
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Sequenced from human survivors of COVID-19 (SARS-CoV-2). The spike NTD is expressed on the surface of SARS-CoV-2.
Preservatives	This monoclonal antibody is aseptically packaged and formulated in 0.01 M phosphate buffered saline (PBS) pH 7.2 - 7.4, 150 mM NaCl with no carrier protein, potassium, calcium or preservatives added.
Form/Appearance	Liquid
Concentration	> 1.0 mg/ml
Storage	Functional grade preclinical antibodies may be stored sterile as received at 2-8° C for up to one month. For longer term storage, aseptically aliquot in working volumes without diluting and store at -80° C. Avoid Repeated Freeze Thaw Cycles.
Note	For research use only
Application notes	ELISA (Quality Tested by Biorbyt); IHC (Other applications as reported in literature)
Isotype	Human IgG1
Clonality	Recombinant
Expiration Date	12 months from date of receipt.



Coating:
Purified
Recombinant
SARS-CoV-
2...