


## Product Datasheet

# SARS-CoV-2 Spike RBD Antibody [2165] (HRP) (orb1274627)

<b>Description</b>	SARS-CoV-2 Spike RBD Antibody [2165] (HRP)
<b>Reactivity</b>	Virus
<b>Conjugation</b>	HRP
<b>Tested Applications</b>	ELISA, IHC
<b>Immunogen</b>	Sequenced from human survivors of COVID-19 (SARS-CoV-2). The spike RBD is expressed on the surface of SARS-CoV-2.
<b>Target</b>	S
<b>Preservatives</b>	This HRP-conjugated antibody is formulated in 0.01 M phosphate buffered saline (PBS) pH 7.4, 1% rAlbumin. (Warning: Use of sodium azide as a preservative will inhibit the enzyme activity of horseradish peroxidase)
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	0.5 mg/ml
<b>Storage</b>	This horseradish peroxidase conjugated monoclonal antibody is stable when stored at 2-8° C. Do not freeze.
<b>Note</b>	For research use only
<b>Application notes</b>	ELISA (Quality Tested by Biorbyt)
<b>Isotype</b>	Human IgG1
<b>Clonality</b>	Recombinant
<b>Expiration Date</b>	12 months from date of receipt.



SARS-CoV-2 Spike RBD Antibody (clone 216...)



Coating:  
Purified  
Recombinant  
SARS-CoV-2...