

Product Datasheet

KMT2D Antibody (orb1277395)

| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody to KMT2D |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IHC |
| Immunogen | Synthetic peptide of human KMT2D |
| Preservatives | pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | Immunogen-specific rabbit IgG |
| Clonality | Polyclonal |
| Uniprot ID | O14686 |
| NCBI | NP_003473 |
| Dilution Range | IHC: 1/50-1/200; ELISA: 1/2000-1/5000 |
| Expiration Date | 12 months from date of receipt. |