

Product Datasheet

B-RAF (phospho-Thr599) antibody (orb127761)

Description Rabbit polyclonal antibody to B-RAF (phospho-Thr599).

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications IHC

Immunogen C-terminal

Concentration 1mg/ml

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.


Note For research use only

Clonality Polyclonal

Uniprot ID **P15056**

Dilution Range IHC-P: 1:50-1:100, ELISA: 1:40000

Expiration Date 12 months from date of receipt.



Immunohistochemical analysis of formalin...