

Product Datasheet

BDKRB2 Antibody (orb1277911)

Description	Rabbit polyclonal antibody to BDKRB2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Synthetic peptide of human BDKRB2
Preservatives	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Immunogen-specific rabbit IgG
Clonality	Polyclonal
MW	44 kDa
Uniprot ID	P30411
NCBI	NP_000614.1
Dilution Range	WB: 1/500-1/2000; IHC: 1/25-1/50; ELISA: 1/5000-1/10000
Expiration Date	12 months from date of receipt.