



Product Datasheet

SIRT3 antibody (orb12785)


Description	Rabbit polyclonal antibody to SIRT3
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthetic peptide of human SIRT3(NP_036371.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	28KDa
Uniprot ID	Q9NTG7
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:20-1:50
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse brain...



Immunohistochemical staining of human li...



Immunohistochemical staining of rat kidn...