



Product Datasheet

IKK-alpha/beta (phospho-Ser176/177) antibody (orb127876)


Description	Rabbit polyclonal antibody to IKK-alpha/beta (phospho-Ser176/177).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Internal
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	O15111, O14920
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from 3...



Western blot analysis of extracts from A...



Immunohistochemistry analysis of formali...