

Product Datasheet

ATP5F1A Antibody (orb1278850)

Description	Rabbit polyclonal antibody to ATP5F1A
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Fusion protein of human ATP5F1A
Preservatives	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Immunogen-specific rabbit IgG
Clonality	Polyclonal
Uniprot ID	P25705
Dilution Range	IHC: 1/50-1/200; ELISA: 1/5000-1/10000
Expiration Date	12 months from date of receipt.