



Product Datasheet

Ephrin B1/B2/B3 (phospho-Tyr324) antibody (orb128003)

Description	Rabbit Polyclonal to Ephrin B1/B2/B3 (phospho-Tyr324).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	C-terminal
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P98172, Q15768, P52799
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:20000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from K...



Immunohistochemical analysis of formalin...