

Product Datasheet

CHMP1A Antibody (orb1280467)

Description	Rabbit polyclonal antibody to CHMP1A
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Fusion protein of human CHMP1A
Preservatives	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Immunogen-specific rabbit IgG
Clonality	Polyclonal
MW	22 kDa
Uniprot ID	Q9HD42
Dilution Range	WB: 1/500-1/2000; IHC: 1/50-1/200; ELISA: 1/2000-1/5000
Expiration Date	12 months from date of receipt.