

## Product Datasheet

### ATP5IF1 Antibody (orb1280915)

<b>Description</b>	Rabbit polyclonal antibody to ATP5IF1
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IHC, WB
<b>Immunogen</b>	Fusion protein of human ATP5IF1
<b>Preservatives</b>	pH 7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	Immunogen-specific rabbit IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	12 kDa
<b>Uniprot ID</b>	<a href="#">Q9UII2</a>
<b>Dilution Range</b>	WB: 1/500-1/2000; IHC: 1/50-1/200; ELISA: 1/5000-1/10000
<b>Expiration Date</b>	12 months from date of receipt.