



Product Datasheet

HYOU1 antibody (orb12812)


Description	Rabbit polyclonal antibody to HYOU1
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human ORP150/GRP170/HSP12A(NP_001124463.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	150kDa
Uniprot ID	Q9Y4L1
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF/ICC: 1:10-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human ve...



Immunohistochemical staining of human co...



Immunohistochemical staining of human co...