

## Product Datasheet

### DLGAP4 Antibody (orb1282238)

<b>Description</b>	Rabbit polyclonal antibody to DLGAP4
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Immunogen</b>	Synthetic peptide of human DLGAP4
<b>Preservatives</b>	pH 7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	Immunogen-specific rabbit IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	108 kDa
<b>Uniprot ID</b>	<a href="#">Q9Y2H0</a>
<b>NCBI</b>	<a href="#">NP_055717</a>
<b>Dilution Range</b>	WB: 1/200-1/1000; ELISA: 1/5000-1/10000
<b>Expiration Date</b>	12 months from date of receipt.