

Product Datasheet

NLRP9 Antibody (orb1283519)

Description	Rabbit polyclonal antibody to NLRP9
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Synthetic peptide of human NLRP9
Preservatives	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Immunogen-specific rabbit IgG
Clonality	Polyclonal
Uniprot ID	Q7RTR0
NCBI	NP_789790
Dilution Range	IHC: 1/50-1/200; ELISA: 1/2000-1/5000
Expiration Date	12 months from date of receipt.