



## Product Datasheet

### MRPL13 antibody (orb128389)

<b>Description</b>	Rabbit Polyclonal to MRPL13.
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC, WB
<b>Immunogen</b>	C-terminal
<b>Concentration</b>	1mg/ml
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<b>Q9BYD1</b>
<b>Dilution Range</b>	WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:40000
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of extracts from A...



Immunohistochemical analysis of formalin...