

Product Datasheet

NSMCE3 Antibody (orb1283974)

Description	Rabbit polyclonal antibody to NSMCE3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Synthetic peptide of human NSMCE3
Preservatives	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Immunogen-specific rabbit IgG
Clonality	Polyclonal
MW	34 kDa
Uniprot ID	Q96MG7
NCBI	NP_619649
Dilution Range	WB: 1/500-1/2000; IHC: 1/50-1/200; ELISA: 1/5000-1/10000
Expiration Date	12 months from date of receipt.