



Product Datasheet

MuSK (phospho-Tyr755) antibody (orb128406)


Description	Rabbit Polyclonal to MuSK (phospho-Tyr755).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	O15146
Dilution Range	WB: 1:500-1:3000, ELISA: 1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of JK cells using ...



Western blot analysis of extracts from J...



Western blot analysis of extracts from J...