

Product Datasheet

SEL1L Antibody (orb1284173)

Description	Rabbit polyclonal antibody to SEL1L
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Synthetic peptide of human SEL1L
Preservatives	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Immunogen-specific rabbit IgG
Clonality	Polyclonal
Uniprot ID	Q9UBV2
NCBI	NP_005056.3
Dilution Range	IHC: 1/25-1/100; ELISA: 1/1000-1/5000
Expiration Date	12 months from date of receipt.