

Product Datasheet

SARS-CoV-2 N Antibody (orb1285891)

Description	Rabbit polyclonal antibody to SARS-CoV-2 N
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Fusion protein corresponding to a region derived from 1-419 amino acids of SARS-CoV-2 N
Preservatives	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Immunogen-specific rabbit IgG
Clonality	Polyclonal
MW	46 kDa
Uniprot ID	PODTC9
Dilution Range	WB: 1/10000-1/50000; ELISA: 1/5000-1/10000
Expiration Date	12 months from date of receipt.