

Product Datasheet

SARS-CoV-2 S2 Antibody (orb1285893)

Description	Rabbit polyclonal antibody to SARS-CoV-2 S2
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Fusion protein corresponding to a region derived from 805-1020 amino acids of SARS-CoV-2 S protein
Preservatives	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Immunogen-specific rabbit IgG
Clonality	Polyclonal
MW	24 kDa
Uniprot ID	P0DTC2
Dilution Range	WB: 1/5000-1/20000; ELISA: 1/5000-1/10000
Expiration Date	12 months from date of receipt.