

## Product Datasheet

### RNF114 Antibody (orb1286414)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Rabbit polyclonal antibody to RNF114  |
| <b>Species/Host</b>        | Rabbit  |
| <b>Reactivity</b>          | Human, Mouse, Rat   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Tested Applications</b> | ELISA, IHC  |
| <b>Immunogen</b>           | Fusion protein of human RNF114  |
| <b>Preservatives</b>       | pH 7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol   |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only   |
| <b>Isotype</b>             | Immunogen-specific rabbit IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Uniprot ID</b>          | <a href="#">Q9Y508</a>  |
| <b>Dilution Range</b>      | IHC: 1/50-1/100; ELISA: 1/5000-1/10000  |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |