



Product Datasheet

eNOS (Phospho-S1177) antibody (orb12872)

Description	Rabbit polyclonal antibody to eNOS (Phospho-S1177)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic phosphorylated peptide around S1177 of human eNOS (NP_000594.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	133kDa
Uniprot ID	P29474
Entrez	4846
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of NIH...



Western blot analysis of extracts of SKO...