



Product Datasheet

STK39 (phospho-Ser311) antibody (orb128793)

Description	Rabbit Polyclonal to STK39 (phospho-Ser311).
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Internal
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9UEW8
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from C...



Immunohistochemical analysis of formalin...