



## Product Datasheet

# Tyrosine Hydroxylase (phospho-Ser19) antibody (orb128835)

<b>Description</b>	Rabbit polyclonal antibody to Tyrosine Hydroxylase (phospho-Ser19).
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IF, IHC, WB
<b>Immunogen</b>	N-terminal
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<b>P07101</b>
<b>Dilution Range</b>	IHC-P: 1:50-1:100, IF/ICC: 1:100-1:500, ELISA: 1:5000
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of extracts from 2...



Immunohistochemistry analysis of formali...