



Product Datasheet

53BP1 (phospho-Ser25) antibody (orb128861)

Description	Rabbit Polyclonal to 53BP1 (phospho-Ser25).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC
Immunogen	N-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q12888
Dilution Range	IHC-P: 1:50-1:100, IF/ICC: 1:100-1:500, ELISA: 1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Confocal immunofluorescence analysis of ...