

Product Datasheet

FLAG Tag Antibody (orb1288617)

Description	Rabbit polyclonal antibody to FLAG Tag
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthetic peptide: DYKDDDDK conjugated to KLH
Preservatives	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Immunogen-specific rabbit IgG
Clonality	Polyclonal
MW	50 kDa
Dilution Range	WB: 1/2000-1/10000; ELISA: 1/5000-1/10000
Expiration Date	12 months from date of receipt.