

Product Datasheet

MAP3K7 (orb1289145)

Description	Mouse monoclonal antibody to MAP3K7
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG2a
Clonality	Monoclonal
Entrez	6885
Expiration Date	12 months from date of receipt.