

## Product Datasheet

### YB1 (phospho-Ser102) antibody (orb128919)

<b>Description</b>	Rabbit Polyclonal to YB1 (phospho-Ser102).	Image not found or type unknown
<b>Species/Host</b>	Rabbit	Western blot analysis of extracts from H...
<b>Reactivity</b>	Human, Mouse, Rat	
<b>Conjugation</b>	Unconjugated	
<b>Tested Applications</b>	WB	
<b>Immunogen</b>	Internal	
<b>Concentration</b>	1mg/ml	
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.	
<b>Note</b>	For research use only	
<b>Clonality</b>	Polyclonal	
<b>Uniprot ID</b>	<b>P67809</b>	
<b>Dilution Range</b>	WB: 1:500-1:3000, ELISA: 1:40000	
<b>Expiration Date</b>	12 months from date of receipt.	