



Product Datasheet

Arginase Microplate Assay Kit (orb1289379)

Description	Arginase Microplate Assay Kit
Reactivity	All
Range	0.1 mmol/L - 3 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Colorimetric
Sample Types	Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	9000-96-8
Expiration Date	Please enquire.



Representative Standard Curve



Positive Control reaction in 96-well pla...