


Product Datasheet

Rat PDGF-BB ELISA Kit (orb1289381)

Description	Rat PDGF-BB ELISA Kit
Reactivity	Rat
Range	62.5 pg/ml - 4000 pg/ml
Tested Applications	ELISA
Target	PDGFB
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q05028
Sensitivity	30 pg/ml
Entrez	24628
Expiration Date	Please enquire.



Representative
Standard
Curve