


Product Datasheet

Mouse Vitamin D-binding protein ELISA Kit (orb1289390)

Description	Mouse Vitamin D-binding protein ELISA Kit
Reactivity	Mouse
Range	312 pg/ml - 20000 pg/ml
Tested Applications	ELISA
Target	GC
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma, Cell Lysates, Tissue Homogenates
Uniprot ID	P21614
Sensitivity	150 pg/ml
Entrez	14473
Expiration Date	Please enquire.



Representative Standard Curve