


Product Datasheet

Mouse Granzyme B ELISA Kit (orb1289415)

Description	Mouse Granzyme B ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml - 1000 pg/ml
Tested Applications	ELISA
Target	GZMB
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P04187
Sensitivity	7 pg/ml
Entrez	14939
Expiration Date	Please enquire.



Representative Standard Curve