


Product Datasheet

Human S100-B ELISA Kit (orb1289441)

| | |
|----------------------------|---|
| Description | Human S100-B ELISA Kit |
| Reactivity | Human |
| Range | 156 pg/ml - 10000 pg/ml |
| Tested Applications | ELISA |
| Target | S100B |
| Storage | Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles. |
| Note | For research use only |
| Application notes | Colorimetric |
| Sample Types | Cell Culture Supernates, Serum, Plasma, Cell Lysates, Tissue Homogenates |
| Uniprot ID | P04271 |
| Sensitivity | 80 pg/ml |
| Entrez | 6285 |
| Expiration Date | Please enquire. |



Representative Standard Curve