


Product Datasheet

Human PIGF ELISA Kit (orb1289446)

| | |
|----------------------------|---|
| Description | Human PIGF ELISA Kit |
| Reactivity | Human |
| Range | 31.2 pg/ml - 2000 pg/ml |
| Tested Applications | ELISA |
| Target | PIGF |
| Storage | Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles. |
| Note | For research use only |
| Application notes | Colorimetric |
| Sample Types | Cell Culture Supernates, Serum, Plasma |
| Uniprot ID | Q07326 |
| Sensitivity | 15 pg/ml |
| Entrez | 5281 |
| Expiration Date | Please enquire. |



Representative
Standard
Curve