


## Product Datasheet

### Human CCL26 ELISA Kit (orb1289488)

|                            |                                                                                               |
|----------------------------|-----------------------------------------------------------------------------------------------|
| <b>Description</b>         | Human CCL26 ELISA Kit                                                                         |
| <b>Reactivity</b>          | Human                                                                                         |
| <b>Range</b>               | 31.2 pg/ml - 2000 pg/ml                                                                       |
| <b>Tested Applications</b> | ELISA                                                                                         |
| <b>Target</b>              | CCL26                                                                                         |
| <b>Storage</b>             | Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles. |
| <b>Note</b>                | For research use only                                                                         |
| <b>Application notes</b>   | Colorimetric                                                                                  |
| <b>Sample Types</b>        | Cell Culture Supernates, Serum, Plasma                                                        |
| <b>Uniprot ID</b>          | <a href="#">Q9Y258</a>                                                                        |
| <b>Sensitivity</b>         | 15 pg/ml                                                                                      |
| <b>Entrez</b>              | <a href="#">10344</a>                                                                         |
| <b>Expiration Date</b>     | Please enquire.                                                                               |



Representative Standard Curve