


Product Datasheet

Human MMP14 ELISA Kit (orb1289496)

Description	Human MMP14 ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10000 pg/ml
Tested Applications	ELISA
Target	MMP14
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P50281
Sensitivity	15 pg/ml
Entrez	4323
Expiration Date	Please enquire.



Representative Standard Curve