

Product Datasheet

Human VEGF165 Protein (orb1289549)

Description	Recombinant Human VEGF165 Protein
Reactivity	Human
Tested Applications	ELISA, MS, SDS-PAGE, WB
Target	VEGFA
Preservatives	Liquid in a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Form/Appearance	Liquid in a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Please minimize freeze-thaw cycles.
Protein Sequence	APMAEGGGQN HHEVVKFMDV YQRSYCHPIE TLVDIFQEYP DEIEYIFKPS CVPLMRCGGC CNDEGLECVP TEESNITMQI MRIKPHQGQH IGEMSFLQHN KCECRPKKDR ARQENPCGPC SERRKHLFVQ DPQTCKCSCK NTDSRCKARQ LELNERTCRC DKPRR
Purity	Greater than 95% as determined by reducing SDS-PAGE.
Source	Human cells
Uniprot ID	P15692
Entrez	7422
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)