


## Product Datasheet

### Ovalbumin antibody (FITC) (orb1289564)

<b>Description</b>	Gallus Polyclonal antibody against Ovalbumin conjugated to FITC
<b>Species/Host</b>	Gallus
<b>Conjugation</b>	FITC
<b>Tested Applications</b>	ELISA, FC, IF, IH
<b>Preservatives</b>	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.02% Sodium Azide
<b>Form/Appearance</b>	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.02% Sodium Azide
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Source</b>	Chicken
<b>Dilution Range</b>	IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)
<b>Expiration Date</b>	12 months from date of receipt.



Line colors represent the approximate vi...