

## Product Datasheet

# Mouse Anti-Chicken IgG (H&L) antibody (Texas Red) (orb1289651)

<b>Description</b>	Mouse Monoclonal antibody against Mouse Anti-Chicken IgG conjugated to Texas Red
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Texas Red
<b>Tested Applications</b>	ELISA, FC, IF, IH
<b>Immunogen</b>	Chicken IgG
<b>Preservatives</b>	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
<b>Form/Appearance</b>	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Monoclonal
<b>Source</b>	Mouse
<b>Dilution Range</b>	IH: 1:100-1:1000, IF: 1:100-1:1000, FC: 1:1000-1:4000, E: Use at an assay dependent concentration
<b>Expiration Date</b>	12 months from date of receipt.