

Product Datasheet

Mouse Anti-Chicken IgG (H&L) antibody (TRITC) (orb1289656)

Description	Mouse Monoclonal antibody against Mouse Anti-Chicken IgG conjugated to TRITC
Species/Host	Mouse
Conjugation	TRITC
Tested Applications	ELISA, FC, IF, IH
Immunogen	Chicken IgG
Preservatives	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
Form/Appearance	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Dilution Range	IH: 1:100-1:1000, IF: 1:100-1:1000, FC: 1:1000-1:4000, E: Use at an assay dependent concentration
Expiration Date	12 months from date of receipt.