



Product Datasheet

COVID-19 S-clgG1 Neutralizing antibody (orb1289751)

Description	Monkey Monoclonal antibody to COVID-19 S-clgG1 Neutralizing
Species/Host	Monkey
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Target	S
Preservatives	1 mg/ml. Cynomolgus IgG1. Liquid in 0.01M Phosphate Buffered Saline, pH 7.4.
Form/Appearance	1 mg/ml. Cynomolgus IgG1. Liquid in 0.01M Phosphate Buffered Saline, pH 7.4.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Cynomolgus
Uniprot ID	P0DTC2
Entrez	43740568
Dilution Range	ELISA: Use at an assay dependent dilution.
Expiration Date	12 months from date of receipt.



Immobilized Recombinant COVID-19 S Prote...



Anti-COVID-19 S-clgG1 Neutralizing Antib...