

Product Datasheet

Hepatitis B Virus Core Antigen-HRP labeled antibody (HRP) (orb1289756)

Description	Goat Polyclonal antibody against Hepatitis B Virus Core Antigen-HRP labeled conjugated to HRP
Species/Host	Goat
Reactivity	Virus
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Hepatitis B Virus Core Antigen
Preservatives	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide.
Form/Appearance	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Goat
Dilution Range	ELISA: Use at an assay dependent dilution.
Expiration Date	12 months from date of receipt.