

Product Datasheet

Hepatitis B Virus Core Antigen antibody (orb1289757)

Description	Goat Polyclonal antibody to Hepatitis B Virus Core Antigen
Species/Host	Goat
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Hepatitis B Virus Core Antigen
Preservatives	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide.
Form/Appearance	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Goat
Dilution Range	ELISA: Use at an assay dependent dilution.
Expiration Date	12 months from date of receipt.