

Product Datasheet

PDGFR beta (phospho-Tyr771) antibody (orb128985)

Description Rabbit polyclonal antibody against PDGFR beta (phospho-Tyr771)

Species/Host Rabbit

Reactivity Human, Rat

Conjugation Unconjugated

Tested Applications WB

Immunogen Internal

Concentration 1mg/ml

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.


Note For research use only

Clonality Polyclonal


Uniprot ID [P09619](#)

Dilution Range WB: 1:500-1:3000, ELISA: 1:1000

Expiration Date 12 months from date of receipt.



Western blot analysis of extracts from H...



Western blot analysis of extracts from c...