

Product Datasheet

Rabbit EGF R Protein (orb1289963)

Description	Rabbit EGF R, His Tag is expressed from human 293 cells (HEK293). It contains AA Leu 25 - Ser 645 (A...
Reactivity	Rabbit
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	EGF R
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 70.2 kDa. The protein migrates as 75-85 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	XP_017199838.1
Purity	95%
Source	Rabbit EGF R, His Tag (orb1289963) is expressed from human 293 cells (HEK293). It contains AA Leu 25 - Ser 645 (Accession # G1SZE4-1).
MW	70.2 kDa
Uniprot ID	G1SZE4
NCBI	XP_017199838.1
Expiration Date	6 months from date of receipt.