



Product Datasheet

P97/DAP5/EIF4G2 Antibody (orb1290016)


Description	P97/DAP5/EIF4G2 Antibody
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IF, WB
Immunogen	E.coli-derived human P97/DAP5/EIF4G2 recombinant protein (Position: E91-E790).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	97 kDa
Uniprot ID	P78344
Expiration Date	12 months from date of receipt.



IF analysis of P97/DAP5/EIF4G2 using ant...



Western blot analysis of P97/DAP5/EIF4G2...



Flow Cytometry analysis of THP-1 cells u...