

## Product Datasheet

### Monkeypox virus A35R antibody (orb1290076)

<b>Description</b>	Mouse monoclonal antibody to Monkeypox virus A35R
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Virus
<b>Endotoxins</b>	Please contact with the lab for this information.
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, NeA
<b>Preservatives</b>	PBS buffer, pH 7.4.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG2a, kappa
<b>Clonality</b>	Monoclonal
<b>Purity</b>	> 95%
<b>Expiration Date</b>	12 months from date of receipt.