

## Product Datasheet

### Monkeypox virus A29L antibody (orb1290077)

<b>Description</b>	Mouse monoclonal antibody to Monkeypox virus A29L
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Virus
<b>Endotoxins</b>	Please contact with the lab for this information.
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, NeA
<b>Preservatives</b>	PBS buffer, pH 7.4.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG2a, lambda
<b>Clonality</b>	Monoclonal
<b>Purity</b>	> 95%
<b>Expiration Date</b>	12 months from date of receipt.