

## Product Datasheet

### Vaccinia L1R/Protein L1 antibody (orb1290085)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Mouse monoclonal antibody to Vaccinia L1R/Protein L1  |
| <b>Species/Host</b>        | Mouse   |
| <b>Reactivity</b>          | Virus   |
| <b>Endotoxins</b>          | Please contact with the lab for this information.   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Tested Applications</b> | ELISA, NeA  |
| <b>Preservatives</b>       | PBS buffer, pH 7.4.   |
| <b>Form/Appearance</b>     | Liquid  |
| <b>Concentration</b>       | 1 mg/ml   |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only   |
| <b>Isotype</b>             | IgG2a, kappa  |
| <b>Clonality</b>           | Monoclonal  |
| <b>Purity</b>              | > 95%   |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |